

Flexible nozzle positioning capability for complete whole room disinfection from inside or outside the room being treated.

The **HaloFogger FLX** provides flexibility in positioning the nozzle while maintaining the precision of disinfectant delivery. The base FLX unit can be located outside the room being treated while the detachable nozzle operates inside the room. The nozzle assembly is detachable and additional hose assemblies can be purchased to facilitate the treatment of multiple spaces quickly and affordably with a single base unit.



Treatment Room Size	up to 9,500 ft ³ / 270 m ³
Timer	10 - 60 minutes
Fluid Output	50 ml (1.7 oz) / minute
Hose Length	12 foot
Target Dose Concentration	11ml/m ³
Average coverage per gallon of disinfectant	11,500 ft ³ / 325 m ³
Average coverage per liter of disinfectant	3,000 ft ³ / 86 m ³



Extra long hose allows the nozzle assembly to be positioned inside a closed room while the base unit remains outside of the sealed door.



Precise dispensing of aerosolized Halo Disinfectant ensures complete coverage of all exposed surfaces with the the room being treated.